



Product Datasheet

Beta-Glucuronidase Microplate Assay Kit (orb707343)

Description	Detection and Quantification of Beta-Glucuronidase Activity.
Reactivity	All
Range	5 µmol/L - 500 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 560 nm
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-45-0
Expiration Date	Please enquire.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...