


Product Datasheet

Sialic Acid Microplate Assay Kit (orb707350)

Description	Detection and Quantification of Sialic Acid Content.
Reactivity	All
Range	10 µmol/L - 1000 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Cell culture media, Saliva, Milk, Tissue extracts, Cell lysate, Other biological fluids, Platform: Microplate reader 570 nm
Sample Types	Serum, Plasma, Cell culture media, Saliva, Milk, Tissue extracts, Cell lysate, Other biological fluids
Expiration Date	Please enquire.



Representative
Standard
Curve