

Product Datasheet

Creatine Microplate Assay Kit (orb707351)

Description Detection and Quantification of Creatine Content.

Reactivity All

Range 0.05 mmol/L - 5 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Other biological fluids, Platform: Microplate reader 546 nm

Sample Types Urine, Serum, Plasma, Other biological fluids

CAS Number 6020-87-7

Expiration Date Please enquire.



Representative
Standard
Curve