


Product Datasheet

Total Carbohydrate Microplate Assay Kit (orb707362)

Description	Detection and Quantification of Total Carbohydrate Content.
Reactivity	All
Range	50 µg/ml - 500 µg/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue extracts, Cell lysate, Food, Juice, Beverage, Other agricultural products, Platform: Microplate reader 540 nm
Sample Types	Serum, Plasma, Tissue extracts, Cell lysate, Food, Juice, Beverage, Other agricultural products
Expiration Date	Please enquire.



Representative
Standard
Curve