

Product Datasheet

Isocitrate Microplate Assay Kit (orb707363)

Description Detection and Quantification of Total Isocitrate Content.

Reactivity All

Range 0.1 mmol/L - 5 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 450 nm

Sample Types Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 320-77-4

Expiration Date Please enquire.



Representative
Standard
Curve