


Product Datasheet

Caspase-9 Microplate Assay Kit (orb707377)

Description	Detection and Quantification of Caspase-9 Activity.
Reactivity	All
Range	0.01 µmol/ml - 5 µmol/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Other biological fluids, Platform: Microplate reader 405 nm
Sample Types	Tissue extracts, Cell lysate, Other biological fluids
Expiration Date	Please enquire.



Representative
Standard
Curve