


Product Datasheet

Xanthine/Hypoxanthine Microplate Assay Kit (orb707387)

Description	Detection and Quantification of Xanthine/Hypoxanthine Content.
Reactivity	All
Range	0.05 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 505 nm
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	69-89-6
Expiration Date	Please enquire.



Representative
Standard
Curve