


Product Datasheet

Ethanol Microplate Assay Kit (orb707393)

Description	Detection and Quantification of Ethanol Content.
Reactivity	All
Range	2 µg/L - 200 µg/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Urine, Saliva, Other biological fluids, Platform: Microplate reader 450 nm
Sample Types	Serum, Plasma, Urine, Saliva, Other biological fluids
CAS Number	64-17-5
Expiration Date	Please enquire.



Representative
Standard
Curve