



Product Datasheet

Glucoamylase Microplate Assay Kit (orb707396)

Description	Detection and Quantification of Glucoamylase Activity.
Reactivity	All
Range	0.1 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Tissue extracts, Cell lysate, Cell culture media, Other biological fluids, Platform: Microplate reader 540 nm
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9032-08-0
Expiration Date	Please enquire.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...