


Product Datasheet

Total Iron-binding Capacity Microplate Assay Kit (orb707397)

Description	Detection and Quantification of Total Iron-binding Capacity.
Reactivity	All
Range	0.1 mmol/L - 10 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Platform: Microplate reader 562 nm
Sample Types	Serum, Plasma
Expiration Date	Please enquire.



Representative
Standard
Curve