

## Product Datasheet

### 2019-nCoV NSP1 protein (orb707413)

|                         |                                                                                                         |
|-------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | Recombinant 2019-nCoV NSP1 is produced by our E.coli expression system and the target gene encoding ... |
| <b>Endotoxins</b>       | Less than 0.1 ng/ug (1 EU/ug)                                                                           |
| <b>Preservatives</b>    | Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.5.              |
| <b>Form/Appearance</b>  | Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.5 .             |
| <b>Storage</b>          | Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.                  |
| <b>Tag</b>              | C-His                                                                                                   |
| <b>Note</b>             | For research use only                                                                                   |
| <b>Protein Sequence</b> | Met1-Gly180                                                                                             |
| <b>Source</b>           | E.coli                                                                                                  |
| <b>Uniprot ID</b>       | YP_009725297.1                                                                                          |
| <b>Expiration Date</b>  | 6 months from date of receipt.                                                                          |