

Product Datasheet

2019-nCov S1 protein protein (orb707439)

Description	Recombinant 2019-nCov S1 protein protein
Preservatives	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-6His
Note	For research use only
Protein Sequence	Val16-Arg685
Source	Human Cells
Uniprot ID	QHD43416.1
Expiration Date	6 months from date of receipt.