

Product Datasheet

SARS-CoV-2 NSP7 protein (orb707445)

Description	Recombinant SARS-CoV-2 NSP7 is produced by our E.coli expression system and the target gene encoding...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 um filtered solution of 20mM Tris, 150mM NaCl, 10%Glycerol, 4mM MgCl2, pH 8.0
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris, 150mM NaCl, 10%Glycerol, 4mM MgCl2, pH 8.0
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-TEV-8His
Note	For research use only
Protein Sequence	Ser1-Gln83
Source	E. coli
Uniprot ID	YP_009725303.1
Expiration Date	6 months from date of receipt.