


Product Datasheet

Tropomyosin (phospho-Tyr674/675) antibody (orb7120)

Description	Tropomyosin (phospho-Tyr674/675) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human TrkA around the phosphorylation site of Tyr674/675 (TD(p-Y)(p-Y)RV)
Target	Phospho-TrkA(Tyr674 + Tyr675) + TrkB(Tyr706 + Tyr707)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90 kDa
Uniprot ID	P04629
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain...