

Product Datasheet

TrkA (phospho-Tyr785) antibody (orb7121)

Description	TrkA (phospho-Tyr785) antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NTRK2 around the phosphorylation site of Tyr817 (PV(p-Y)LD)
Target	Phospho-TrkB (Tyr817)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90 kDa
Uniprot ID	Q16620
Expiration Date	12 months from date of receipt.