

## Product Datasheet

### Human Activated Protein C (390-404) peptide (orb71342)

|                         |                                                             |
|-------------------------|-------------------------------------------------------------|
| <b>Description</b>      | Human Activated Protein C (390-404) peptide                 |
| <b>Form/Appearance</b>  | White lyophilized powder                                    |
| <b>Storage</b>          | Keep in cool and dry place                                  |
| <b>Note</b>             | For research use only                                       |
| <b>Protein Sequence</b> | Tyr-Gly-Val-Tyr-Thr-Lys-Val-Ser-Arg-Tyr-Leu-Asp-Trp-Ile-His |
| <b>Purity</b>           | > 95%                                                       |
| <b>Expiration Date</b>  | 12 months from date of receipt.                             |