

Product Datasheet

Endothelin 3 peptide (orb71562)

Description	Endothelin 3 peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Cys-Thr-Cys-Phe-Thr-Tyr-Lys-Asp-Lys-Glu-Cys-Val-Tyr-Tyr-Cys-His-Leu-Asp-Ile-Ile-Trp(Disulfide bridge Cys1-Cys15,Cys3-Cys11)
Purity	> 95%
Expiration Date	12 months from date of receipt.