

Product Datasheet

Histone H3 (116136), C116136 peptide (orb71672)

Description	Histone H3 (116136), C116136 peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Lys-Arg-Val-Thr-Ile-Met-Pro-Lys-Asp-Ile- Gln-Leu-Ala-Arg-Arg-Ile-Arg-Gly-Glu- Arg-Ala
Purity	> 95%
Expiration Date	12 months from date of receipt.