

Product Datasheet

Histone H3 (116136), N1539 peptide (orb71673)

Description	Histone H3 (116136), N1539 peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Ala-Pro-Arg-Lys-Gln-Leu-Ala-Thr-Lys-Ala-Ala-Arg-Lys-Ser-Ala-Pro-Ala-Thr-Gly-Gly-Val-Lys-Lys-Pro-His
Purity	> 95%
Expiration Date	12 months from date of receipt.