

Product Datasheet

Human Hypocretin (70-98) peptide (orb71705)

Description	Human Hypocretin (70-98) peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Arg-Ser-Gly-Pro-Pro-Gly-Leu-Gln-Gly-Arg- Leu-Gln-Arg-Leu-Leu-Gln-Ala-Ser-Gly- Asn-His-Ala-Ala-Gly-Ile-Leu-Thr-Met-Gly
Purity	> 95%
Expiration Date	12 months from date of receipt.