

Product Datasheet

Human I-309 peptide (orb71706)

Description	Human I-309 peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Lys-Ser-Met-Gln-Val-Pro-Phe-Ser-Arg-Cys-Cys-Phe-Ser-Phe-Ala-Glu-Gln-Glu-Ile-Pro-Leu-Arg-Ala-Ile-Leu-Cys-Tyr-Arg-Asn-Thr-Ser-Ser-Ile-Cys-Ser-Asn-Glu-Gly-Leu-Ile-Phe-Lys-Leu-Lys-Arg-Gly-Lys-Glu-Ala-Cys-Ala-Leu-Asp-Thr-Val-Gly-Trp-Val-Gln-Arg-His-Arg-Lys-Met
Purity	> 95%
Expiration Date	12 months from date of receipt.