

Product Datasheet

LL37 pent peptide (orb71736)

Description	LL37 pent peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Leu-Leu-Gly-Asn-Phe-Phe-Arg-Lys-Ser-Lys-Gln-Lys-Ile-Gly-Lys-Gln-Phe-Lys-Arg-Ile-Val-Gln-Arg-Ile-Lys-Asn-Phe-Phe-Arg-Asn-Leu-Val-Pro-Arg-Thr-Gln-Ser
Purity	> 95%
Expiration Date	12 months from date of receipt.