

Product Datasheet

Pig Neuromedin (U25) peptide (orb71790)

Description	Pig Neuromedin (U25) peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Phe-Lys-Val-Asp-Glu-Glu-Phe-Gln-Gly-Pro-Ile-Val-Ser-Gln-Asn-Arg-Arg-Tyr-Phe-Leu-Phe-Arg-Pro-Arg-Asn-NH2
Purity	> 95%
Expiration Date	12 months from date of receipt.