

Product Datasheet

Rat Neuromedin S peptide (orb71792)

Description	Rat Neuromedin S peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Leu-Pro-Arg-Leu-Leu-His-Thr-Asp-Ser-Arg-Met-Ala-Thr-Ile-Asp-Phe-Pro-Lys-Lys-Asp-Pro-Thr-Thr-Ser-Leu-Gly-Arg-Pro-Phe-Phe-Leu-Phe-Arg-Pro-Arg-Asn-NH ₂
Purity	> 95%
Expiration Date	12 months from date of receipt.