

Product Datasheet

Bovine Pancreatic Poly peptide (orb71847)

Description	Bovine Pancreatic Poly peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Ala-Pro-Leu-Glu-Pro-Glu-Tyr-Pro-Gly-Asp- Asn-Ala-Thr-Pro-Glu-Gln-Met-Ala-Gln-Tyr- Ala-Ala-Glu-Leu-Arg-Arg-Tyr-Ile-Asn-Met- Leu-Thr-Arg-Pro-Arg-Tyr-NH ₂
Purity	> 95%
Expiration Date	12 months from date of receipt.