

Product Datasheet

Human Pahyroid Hormone (1-34)-Lys peptide (orb71855)

Description	Human Pahyroid Hormone (1-34)-Lys peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Ser-Val-Ser-Glu-Ile-Gln-Leu-Met-His-Asn-Leu-Gly-Lys-His-Leu-Asn-Ser-Met-Glu-Arg-Val-Glu-Trp-Leu-Arg-Lys-Lys-Leu-Gln-Asp-Val-His-Asn-Phe-Lys(Biotin)
Purity	> 95%
Expiration Date	12 months from date of receipt.