

## Product Datasheet

# Phosphorylated Protein Kinase C Substrate 2 peptide (orb71880)

|                         |                                                     |
|-------------------------|-----------------------------------------------------|
| <b>Description</b>      | Phosphorylated Protein Kinase C Substrate 2 peptide |
| <b>Form/Appearance</b>  | White lyophilized powder                            |
| <b>Storage</b>          | Keep in cool and dry place                          |
| <b>Note</b>             | For research use only                               |
| <b>Protein Sequence</b> | Lys-Arg-pThr-Ile-Arg-Arg                            |
| <b>Purity</b>           | > 95%                                               |
| <b>Expiration Date</b>  | 12 months from date of receipt.                     |