

Product Datasheet

Protein Kinase C (530-558) peptide (orb71929)

Description	Protein Kinase C (530-558) peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Leu-Leu-Tyr-Glu-Met-Leu-Ala-Gly-Gln-Ala-Pro-Phe-Glu-Gly-Glu-Asp-Glu-Asp-Glu-Leu-Phe-Gln-Ser-Ile-Met-Glu-His-Asn-Val-NH ₂
Purity	> 95%
Expiration Date	12 months from date of receipt.