

Product Datasheet

Human pTH-Related Protein Splice Isoform 3 (140-173) peptide (orb71933)

Description	Human pTH-Related Protein Splice Isoform 3 (140-173) peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Thr-Ala-Leu-Leu-Trp-Gly-Leu-Lys- Lys-Lys-Lys-Glu-Asn-Asn-Arg-Arg- Thr-His-His-Met-Gln-Leu-Met-Ile- Ser-Leu-Phe-Lys-Ser-Pro-Leu- Leu-Leu-Leu
Purity	> 95%
Expiration Date	12 months from date of receipt.