

Product Datasheet

TAT 2-4 peptide (orb71993)

Description	TAT 2-4 peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Tyr-Gly-Arg-Lys-Lys-Arg-Arg-Gln-Arg-Arg-Arg-Gly-Tyr-Gly-Arg-Lys-Lys-Arg-Arg-Gln-Arg-Arg-Arg-Gly
Purity	> 95%
Expiration Date	12 months from date of receipt.