

Product Datasheet

Urocortin 3 peptide (orb72021)

Description	Urocortin 3 peptide
Form/Appearance	White lyophilized powder
Storage	Keep in cool and dry place
Note	For research use only
Protein Sequence	Phe-Thr-Leu-Ser-Leu-Asp-Val-Pro-Thr-Asn-Ile-Met-Asn-Ile-Leu-Phe-Asn-Ile-Asp-Lys-Ala-Lys-Asn-Leu-Arg-Ala-Lys-Ala-Ala-Ala-Asn-Ala-Gln-Leu-Met-Ala-Gln-Ile-NH ₂
Purity	> 95%
Expiration Date	12 months from date of receipt.