


Product Datasheet

CCP peptide (orb72449)

Description	CCP peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	CCP
Form/Appearance	Lyophilised
Storage	Store at -20 °C, avoid repeated freeze/thaw cycles.
Note	For research use only
Protein Sequence	HQCHQEST citrulline GRSRGRCGRSGS; disulfide bonds: C3 = C16
Source	Synthetic
Expiration Date	6 months from date of receipt.



Auto-Scaled
Chromatogram



CCP peptide