

Product Datasheet

PRKAA2 (phospho-Ser173) antibody (FITC) (orb7261)

Description	PRKAA2 (phospho-Ser173) antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human AMPK alpha 2 around the phosphorylation site of Ser173
Target	phospho-AMPK alpha 2 (Ser173)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P54646
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.