

Product Datasheet

Insulin Receptor alpha peptide [Out of stock] (orb72875)

Description	Insulin Receptor alpha peptide.
Target	INSR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb10923 (Insulin Receptor alpha antibody)
CAS Number	3
Uniprot ID	P06213
Expiration Date	12 months from date of receipt.