



Product Datasheet

ZNF398 antibody (orb73389)

Description	ZNF398 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q8TD17
Expiration Date	12 months from date of receipt.



Western Blot
analysis of
A:293T



B:H1299
whole cell
lysate using
ZNF398 a...