

Product Datasheet

RV-C 1-123 peptide (orb73400)

Description	The protein contains the Rubella Virus Capsid C region, amino acids 1-123. The protein is purified b...
Form/Appearance	liquid (Supplied in 25 mM Tris-HCl pH 8.0, 10 mM Glutathione, 0.2% Sakosil, 50% Glycerol); pH: 8.0
Storage	store at -20 °C
Note	For research use only
Application notes	Antigen in ELISA and Western blots, excellent antigen for detection of Rubella virus with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE, RP-HPLC)
Expiration Date	12 months from date of receipt.