



Product Datasheet

SERPINI2 antibody (orb73407)

Description	SERPINI2 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	075830
Expiration Date	12 months from date of receipt.



Western Blot analysis of A:White adipose...



Immunohistochemical analysis of formalin...