

Product Datasheet

ABCF2 Polyclonal Antibody [Out of stock] (orb734246)

Description	ABCF2 Polyclonal Antibody
Species/Host	Rabbit
Tested Applications	IF, IHC-P, WB
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.190.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UG63
Expiration Date	12 months from date of receipt.