

Product Datasheet

HIV-1 Nef 3-190 peptide (orb73426)

Description	The protein contains the fragment of HIV-1 Nef protein, amino acids 3-190, the molecular weight is 2...
Form/Appearance	liquid (Supplied in 20 mM Tris-HCl pH 8.0, 8 M urea and 10 mM beta-mercaptoethanol)
Storage	store at -20 °C
Note	For research use only
Application notes	May be used in ELISA and Western blots, excellent antigen for early detection of HIV seroconvertors with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE, HPLC)
MW	Theoretical MW: 134 kDa
Expiration Date	12 months from date of receipt.