

Product Datasheet

TG MIC3 protein (orb73432)

Description	The protein contains the Toxoplasma gondii MIC3 immunodominant regions, amino acids 234-306. The pro...
Form/Appearance	liquid (Supplied in 50 mM Tris-HCl pH 8.0, 1.5 M urea and 0, 1 mM EDTA)
Storage	store at -20 °C
Note	For research use only
Application notes	Antigen in ELISA and Western blots, excellent antigen for detection of Toxoplasma gondii with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE)
Source	E. coli
Expiration Date	12 months from date of receipt.