

Product Datasheet

CMV pp28 protein (orb73434)

Description	The protein contains the CMV pp28 immunodominant regions, amino acids 130-160. The protein is purifi...
Form/Appearance	liquid (Supplied in 50 mM Tris-HCl pH 7.2, 1 mM EDTA and 50% glycerol)
Storage	store at -20 °C
Note	For research use only
Application notes	Antigen in ELISA and Western blots, excellent antigen for detection of CMV with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE, RP-HPLC)
Expiration Date	12 months from date of receipt.