



Product Datasheet

Cortisol Binding Globulin antibody (orb73445)

Description	Cortisol Binding Globulin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P08185
Expiration Date	12 months from date of receipt.



Western Blot
analysis of A:A431
whole ce...



Immunofluorescence
analysis of formalin-

...