

## Product Datasheet

### Toxoplasma gondii GRA1/p24 57-149 protein (orb73458)

|                          |                                                                                                                             |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | The protein contains the p24 (GRA1) immunodominant regions, amino acids 57-149. The protein is purif...                     |
| <b>Form/Appearance</b>   | liquid (Supplied in 1x PBS pH 7.2 and 50% glycerol)                                                                         |
| <b>Storage</b>           | store at -20 °C                                                                                                             |
| <b>Note</b>              | For research use only                                                                                                       |
| <b>Application notes</b> | Antigen in ELISA and Western blots, excellent antigen for detection of Toxoplasma gondii with minimal specificity problems. |
| <b>Clonality</b>         | Recombinant                                                                                                                 |
| <b>Purity</b>            | > 90% (SDS-PAGE)                                                                                                            |
| <b>Source</b>            | E. coli                                                                                                                     |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                             |