

Product Datasheet

Mouse CLEC12A / M1CL / CLL-1 Protein (orb735038)

Description	Mouse CLEC12A, His Tag is expressed from human 293 cells (HEK293). It contains AA Tyr 65 - Arg 267 (...)
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CLEC12A
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 25.7 kDa. The protein migrates as 30-40 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Tyr 65 - Arg 267
Purity	90%
Source	Mouse CLEC12A, His Tag (orb735038) is expressed from human 293 cells (HEK293). It contains AA Tyr 65 - Arg 267 (Accession # Q504P2-1).
MW	25.7 kDa
Uniprot ID	Q504P2
Expiration Date	6 months from date of receipt.