

Product Datasheet

Biotinylated Human IL-22 Protein (orb735043)

Description

Biotinylated Human IL-22, His,Avitag (IL2-H8247) is

Biotinylated Human IL-22, His,Avitag (IL2-H8247) is

Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Biotin
Target	IL-22
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	N-10×His & N-Avi
Note	For research use only

Application notes This protein carries a polyhistidine tag at the N-terminus, followed by an Avi tag (). The protein has a calculated MW of 20.3 kDa. The protein migrates as 25-35 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity


- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source Biotinylated Human IL-22, His, (orb735043) is expressed from human 293 cells (HEK293). It contains AA Ala 34 - Ile 179 (Accession # Q9GZX6-1).

MW 20.3 kDa

Uniprot ID [Q9GZX6](https://www.uniprot.org/entry/Q9GZX6)

Expiration Date 6 months from date of receipt.



Immobilized
Biotinylated
Human IL-
22, Hi...