

Product Datasheet

Human MIS RII Protein (orb735045)

Description	MIS RII (Mullerian inhibiting substance type II receptor), also known as AMHRII (anti-Mullerian horm...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	MIS RII
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 40.0 kDa. The protein migrates as 50-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Human MIS RII, Fc Tag (orb735045) is expressed from human 293 cells (HEK293). It contains AA Pro 18 - Ser 144 (Accession # Q16671-1).
MW	40.0 kDa
NCBI	NP_065434.1
Expiration Date	6 months from date of receipt.