

## Product Datasheet

### Human IL24 ELISA Kit (orb735144)

|                            |                       |
|----------------------------|-----------------------|
| <b>Description</b>         | Human IL24 ELISA Kit  |
| <b>Reactivity</b>          | Human                 |
| <b>Range</b>               | 62.50 - 4,000 pg/mL   |
| <b>Tested Applications</b> | ELISA                 |
| <b>Note</b>                | For research use only |
| <b>Sensitivity</b>         | 31 pg/mL              |
| <b>Expiration Date</b>     | Please enquire.       |