

## Product Datasheet

### Human Park7 ELISA Kit (orb735185)

|                            |                       |
|----------------------------|-----------------------|
| <b>Description</b>         | Human Park7 ELISA Kit |
| <b>Reactivity</b>          | Human                 |
| <b>Range</b>               | 62.50 - 4000 pg/mL    |
| <b>Tested Applications</b> | ELISA                 |
| <b>Note</b>                | For research use only |
| <b>Sensitivity</b>         | 32 pg/mL              |
| <b>Expiration Date</b>     | Please enquire.       |