

Product Datasheet

RV gE1 157-176/374-390/213-239 peptide (orb73531)

Description	The protein contains glycoprotein E1 immunodominant fragments, amino acids: 157-176, 213-239, and 37...
Form/Appearance	liquid (Supplied in 1x PBS)
Storage	store at -20 °C
Note	For research use only
Application notes	Antigen in ELISA and Western blots, excellent antigen for detection of Rubella virus with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE, RP-HPLC)
Expiration Date	12 months from date of receipt.