


Product Datasheet

TGF beta induced factor 2 antibody (orb73557)

Description	TGF beta induced factor 2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9GZN2
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
A:A431,
B:Hela ...